













19-20 Berners Stree	et, London W1T 3NW	020 7387 6930 E	mail: contact@nimr	odental.co.uk
Prescription No:	Date:		E-mail:	
Clinician:	Patient's Name:	· -		
Postcode: Clinician's A		ress:		
Date Required:				
■ EMERGENCY SERVICE FOR APPLIANCES* ■ PATIENT WILL BE APPOINTED AFTER APPLIANCE ARRIVES ■ RETURN DUPLICATE SET OF MODE APPLIANCE INSURANCE		APPLIANCES*  Powerprox Splint (includes 3x3 and Essix)  3x3 (Fixed Bonded Retainer)  Essix Inman Aligners Spring Retainer  Hawley Myofunctional Tongue Bead Nitanium Molar Rotator Nitanium Palatal Expander FEES APPLY		
INDIRECT BONDING SYSTEINCLUdes: Digital Study Models with IPF (See Below for Details), plus Indirect Bonding Tray(s) Composite Upper AND Lower Upper Upper UD 018 Bracket UR2 Bracket for Shallow/Open Bite INDIRECT BONDING SYSTEIN Indirect Bonding Tray(s) Composite Upper AND Lower Upper Digital Study Models with IPR Analyst Individual Patient Kit (See Below for I	R Analysis and Individual Patient Kit Brackets Stainless Steel Brackets Lower UL2  M - INDIVIDUAL REQUEST OPTIONS Brackets Stainless Steel Brackets Lower sis	UPPER   FIXED   REMOVABLE	- COCOCA (	
Custom Height Standard Hoight  3.5mm	4.0mm 4.5mm 5.0mm 4.0mm 4.5mm	4.5mm 4.0mm 5.0mm 4.0mm 4.0mm 5.0mm 4.5mm	n 4.0mm 3.5mm F	Custom Height Standard Height L L Indiandard Height Heigh
	IATURE			
GO GREEN! PLEASE SCAN OR MAKE A COPY OF THIS PRESCRIPTION FORM FOR YOUR RECORDS  Individual Patient Kit  Brackets: Composite Brackets (1 set) - Upper and Lower				

Archwires: One each of the following

Upper .014 Tooth-Colored NiTi Wire Lower .014 Tooth-Colored NiTi Wire Upper .016 Tooth-Colored NiTi Wire Lower .016 Tooth-Colored NiTi Wire Upper .018 Tooth-Colored NiTi Wire Lower .018 Tooth-Colored NiTi Wire

Upper .018 Tooth-Colored Stainless Steel Wires Lower .018 Tooth-Colored Stainless Steel Wires Upper 16x22 Tooth-Colored Stainless Steel Wires Lower 16x22 Tooth-Colored Stainless Steel Wires

Upper 016 NiTi Reverse Curve - Straight Leg

Auxillaries: Clear Ligature Ties (1 pack), Clear Chain Elastics (1 roll), Bracket Bonding Agent Ortho Solo (1 set-up), Bracket Adhesive Enlight (1 set-up), Diamond Disc w/ Clear Guard (1), Diamond Strip (1), Satin Diamond Bur (1), Shorty

Twist Ties (pk 10), Patient Wax (1 each).